



Wall-mounted

Tabletop

- Crystal clear medical image
- Fast-upto 70 full-size plates per hour
- Ultra light & compact
- Strong and robust cassette
- Accommodates 3 standard medical cassette
- One step, integrated scanning & erasing
- Quiet, stable & low-maintenance



IP cassette

Fire CR Flash Medical Reader

Specification

Sampling Pixel Pitch	Standard	200um
	High	100um
Pixel Matrix (IP 35cm x 43cm)	Standard	1750 x 2150
	High	3500 x 4300
Pixel Matrix (IP 24cm x 30cm)	Standard	1200 x 1500
	High	2400 x 3000
Pixel Matrix (IP 18cm x 24cm)	Standard	900 x 1200
	High	1800 x 2400
Accepted Cassette Size	35cm x 43cm, 24cm x 30cm, 18cm x 24cm	
AD Convert Resolution	16 bit	
Eraser	Embedded	
Computer Interface	USB 2.0 / Ethernet	
Dimensions	456 x 803 x 137mm	
	17.9 x 31.6 x 5.4inch	
Weight	19.5kg	
	42.4lbs	
Power Requirement	100 ~ 240V / 50 ~ 60Hz	
System Configuration	Tabletop or Wall Mount	
Image File Format	DICOM 3.0	

Image processing S/W

QuantorMed Plus –




New Study
Patient information,
Add new study



Region selection
Select the region to
scan



Scan
Insert a cassette and
scan the image



Review
Useful conducting
diagnosis the image